



Bring professional & enjoyable ear cleaning home

Open bebird D3 Pro

IOS
Search for "bebird" in App Store
Android
Search for "bebird" in Google Play



1

Start bebird D3 Pro
Scan QR code to
download APP




2

Open APP
Click bebird D3 Pro

3

Pull to turn on
automatically



Connect Product WiFi


4

Connect Product Wi-Fi
"bebird-D3 Pro-XXXX"



5

Back to bebird APP
Enter and start to use it



Accessory Replacement

Spinning to open the fixing ring to replace the ear spoon. Insert the buckle of ear spoon, spinning to tighten the fixing ring

Insert and fasten the ear spoon head



The ear spoon must be assembled in place before using to prevent falling off

Charging slot



Bottom Storage Box

Spin to open the charging base to swap accessories



Accessory

1

Metal & Silicone wrapped Ear Spoon (Small) x3
Harmless ear cleaning



2

Microlens Cover x 1
Cover the magnifying target with the Plug the ear cleaning rod in the lens cover, it turns into an electronic magnifier



2

Double-color Q-elastic Ear Spoon (Small) x 2
Thorough removing of cerumen



6

Sponge Ear Spoon x 3
Removing of earwax and grease in ears



3

Double-color Annular Ear Spoon x 2
Easy Scraping; higher curvature



7

Charging Cable x 1
TYPE-C Charging cable




4

Double-color Q-elastic Ear Spoon (Big) x 2
Easy and thorough removing of cerumen



Product Size



Hazardous substance

Part name	Hazardous material					
	(Pb)	(Hg)	(Cd)	(Cr(VI))	(PBB)	(PBDE)
Printed Circuit Components	×	○	○	○	○	○
Case	○	○	○	○	○	○
Plastic rubber accessory	○	○	○	○	○	○

This form is formulated in accordance with SJ/T 11364
 ○ : represent the hazardous material ratio is below GB/T26572-2011
 × : represent the hazardous material ratio is over GB/T26572-2011

Troubleshooting

Trouble	Solution
starting up and no response	Please use the attached charging cable for quick charging. The charging time is 1hr
Camera does not turn on	Close and open the cover to try a few times
APP quit quickly	Delete the existing APP and re-scan the QR code for installation
Installation of tools in accessory barrel is not smooth	Firstly align the ear-cleaning tool with the lens sleeve, then install it into the lens sleeve
The power indicator of the ear picker shows blue light Flashes quickly	Please use the attached charging cable for charging. The charging time is 1hr
The display on the app is not reading properly	Disconnect the phone with product and open cover to restart it, then set phone to reconnect product
WiFi connection is unstable, sometimes it is disconnected in seconds	Close and open the cover to restart rod and restart app in phone to connect rod

If you still have problems with the above procedures, please contact our Customer Service Representatives directly for further assistance

Attention

- When cleaning lens, please wipe it carefully with professional alcohol cotton swab.
- Please don't use it while you are moving, and protect from being hit by others.
- This product is unsuitable for children under 3 year old.
- It is forbidden for children to use this product alone in order to avoid accidental harm.
- Please do not put the product in liquid to avoid damage.
- Product built-in rechargeable lithium battery, so it must be charged once a month if it is not used frequently.
- Avoid product exposure to sunshine directly, especially the accessories, in order to avoid softening.
- During using it, the temperature of the product will rise slightly (up to 35 degree), please rest assured to use.
- To ensure the best experience, please keep the child in a sitting position and operate the product horizontally.



Product parameters

Brand	bebird
Product Name	Smart visual ear-clean Rod
Model	D3 Pro
Rod weight	120g
Network standard	IEEE 802.11b/g/n
Antenna	Internal FPC antenna
Working frequency	2.4Ghz
Image transmission rate	20fps
image sensor	CMOS
Working temperature	-10 ~ 45°C
Battery	230mAh lithium battery
Battery life	about 70 minutes
Charging time	1h
Input Parameters	5V ~ 0.4A
Working Voltage	4.2V
Amount of Power	5W
Lens diameter	3.5mm
Optimum focal length	1.5 ~ 2cm
Depth of field	10 ~ 50mm
Pixel	8M
Gyroscope	360 degree multi-directional gyroscope



Warranty card

User information	
User name	Telephone
Shipping Address	
Sales company	Telephone



bebird

Bring professional & enjoyable ear cleaning home

FCC Warning

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Radiation Exposure Statement

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.